

I hereby certify that this correspondence is being deposited with the United States Postal Service as first class mail in an envelope addressed to:

Commissioner for Patents
P. O. Box 1450
Alexandria, VA 22313-1450

On December 17, 2003
TOWNSEND AND TOWNSEND AND CREW
LLP

By: *Dana Kane*
Dana Kane

PATENT

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re application of:

Michael I. Watkins and Richard B.
Edwards

Application No.: 09/905,338

Filed: July 13, 2001

For: MULTIPLEX FLOW ASSAYS
PREFERABLY WITH MAGNETIC
PARTICLES AS SOLID PHASE

Examiner: Stucker

Art Unit: 1648

SECOND DECLARATION UNDER 37
C.F.R. 1.131

Assistant Commissioner for Patents
Washington, D.C. 20231

MICHAEL I. WATKINS and RICHARD B. EDWARDS declare and state:

1. We are the inventors of the invention claimed in claims 21-29 and 50-58 of this Application.
2. The attached exhibit A is a photocopy of laboratory notebook entries and other materials describing experimental work that was carried out in the United States, a NAFTA country or a WTO country.
3. The experimental work described in Exhibit A was conducted prior to March 14, 1997.

4. The experimental work described in the attached Exhibit A was carried out by one or both of us, or by a person acting under the supervision of one or both of us.
5. The experimental work described in the attached Exhibit A corresponds to Examples 1 and 3 of this Application, and shows an experiment in which a plurality of types of magnetic beads was used to detect multiple analytes in a sample using flow cytometry.
6. As shown in Exhibit A, three types of beads were utilized - two sizes of SPHERO™ Carboxyl magnetic particles and one type of SINTEF™ magnetic particles. The three types of beads were differentiable from one another by particle size subrange. Each group of beads was combined with a different antigen.
7. As shown in Exhibit A and described in this patent application, the three types of beads were:

SPHERO™ Carboxyl Magnetic particles, from Spherotech, Inc.,

Libertyville, Illinois, USA -- poly(styrene/acrylic acid particles),
4.35 micrometers (μm) in diameter, density 1.17 g/cc, containing
12% magnetite (by weight)

SPHERO™ Carboxyl Magnetic particles, from Spherotech, Inc.,

Libertyville, Illinois, USA -- poly(styrene/acrylic acid particles),
3.18 μm in diameter, density 1.17 g/cc, containing 12% magnetite
(by weight)

SINTEF Applied Chemistry, Trondheim, Norway --

poly(styrene/divinylbenzene) particles, 10 μm in diameter, density
1.23 g/cc, containing 17.9% magnetite/maghemite (by weight)

8. As shown in Exhibit A, pp. 1, 3 and 5, the particles were coupled to CMV, HSV2 and RUB antigens, respectively. Pages 2, 4 and 6 describe the beads. As shown in Exhibit A, pp. 7 and 8, the particles were then mixed and contacted with patient samples having known quantities of CMV, HSV2 and RUB antigens, including combinations of such

antigens, and were subjected to flow cytometry. The results are shown in the table on p. 7 of Exhibit A and below in Table II, demonstrating that multiple analytes could be detected using the magnetic particles described, in a flow cytometric immunoassay. Page 8 of Exhibit A

9. More specifically, the experimental procedure shown in the attached Exhibit A was as follows:

TABLE I
Amounts Used

Bead	Viral Antigen	Amount of Beads	Weight of Viral Antigen	Volume of Viral Antigen	Volume of Phosphate Buffer (100 mM)
4.35 μm	CMV	10 mg	225.8 μg	322.6 μL	677.4 μL
3.18 μm	HSV2	5 mg	163.0 μg	815.0 μL	185.0 μL
10 μm	RUB	5 mg	5.2 μg	104.0 μL	896.0 μL

The beads in each case were placed in test tubes and washed multiple times with 100 mM phosphate buffer, pH 6.8. The washed beads were then suspended in the volume of phosphate buffer listed in Table I, and respective antigen solution was added (CMV antigen from Chemicon International Incorporated, Temecula, California, USA; HSV2 antigen from Ross Southern Labs, Salt Lake City, Utah, USA; and RUB antigen from Viral Antigens, Memphis, Tennessee, USA) in the amount listed in Table 1. The test tubes were then rotated in end-over-end fashion overnight at room temperature. The tubes were then placed on a magnetic separator and the supernatant was drawn off and discarded. The resulting beads were washed with a wash buffer consisting of 50 mM phosphate buffer, pH 7.4, 0.01% Tween 20, 1% bovine serum albumin, 0.1% sodium azide, 150 mM sodium chloride, then again subjected to magnetic separation, and suspended in a storage buffer consisting of 50 mM phosphate buffer, pH 7.4, 5% glycerol, 1% bovine serum albumin, 0.1% sodium azide, 150 mM sodium chloride.

Procedure:

1. 100 μ L each of five of patient samples (diluted 1:10 in wash buffer), of known CMV, HSV2 and RUB antibody status, were added to 12 \times 75 mm polypropylene test tubes.
2. To each tube was added 100 μ L of a mixture of CMV, HSV2 and RUB antigen-coated particles (described in Example 1) diluted in wash buffer.
3. The tubes were vortexed at ambient temperature for 15 minutes.
4. After vortexing, 800 μ L of wash buffer was added to each tube.
5. The tubes were placed in a magnetic separator for 5 minutes and the liquid phase removed.
6. Steps 4 and 5 were repeated, but with 1000 μ L of wash buffer.
7. 200 μ L of a 1:300 dilution of anti-human IgG-phycoerythrin conjugate (Chemicon International Inc., Temecula, California, USA) was added.
8. The tubes were vortexed at ambient temperature for 15 minutes.
9. After this time, the samples were injected into a flow cytometer (Bryte HS, Bio-Rad Laboratories, Inc., Hercules, California, USA) equipped with a xenon arc lamp.

The results are summarized in Table II below. The data show that the positive samples had increased fluorescence relative to the negative samples. Testing of samples containing only RUB shows that essentially the same results are obtained for a particular sample whether it is assayed with only one particle size directed towards a single analyte (RUB) or with particles of different sizes, each size being directed towards a different analyte.

TABLE II

Test Results

Sample	Antibody Status			Relative Linear Fluorescence Units		
	CMV	HSV2	RUB	CMV	HSV2	RUB
CN6	+	-	+	14	7	155
CN8	+	-	+	16	6	181
CN12	-	-	+	5	7	240
CN15	-	-	+	5	6	329
23	-	+	-	5	45	43

We further declare that all statements made herein of our own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under 18 U.S.C. § 1001, and that such willful false statements may jeopardize the validity of the patent to which this verified states is directed.

Michael I. Watkins

Date: _____



Richard B. Edwards

Date: Dec 1, 2003

TOWNSEND and TOWNSEND and CREW LLP
Two Embarcadero Center, 8th Floor
San Francisco, California 94111-3834
Tel: (415) 576-0200
Fax: (415) 576-0300
JA:ja

TABLE II

Test Results

Sample	Antibody Status			Relative Linear Fluorescence Units		
	CMV	HSV2	RUB	CMV	HSV2	RUB
CN6	+	-	+	14	7	155
CN8	+	-	+	16	6	181
CN12	-	-	+	5	7	240
CN15	-	-	+	5	6	329
23	-	+	-	5	45	43

We further declare that all statements made herein of our own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under 18 U.S.C. § 1001, and that such willful false statements may jeopardize the validity of the patent to which this verified states is directed.


Michael I. Watkins

Date: 10/20/03

Richard B. Edwards

Date: _____

TOWNSEND and TOWNSEND and CREW LLP
Two Embarcadero Center, 8th Floor
San Francisco, California 94111-3834
Tel: (415) 576-0200
Fax: (415) 576-0300
JA:ja

Adsorption of CMV Antigen to Magnetic Beads

Purpose: To adsorb ~~differ~~ Chemicon CMV antigen to magnetic beads.

Procedure

Tube	Bead	Amt Beads	Vol. Beads	Vol. CMV ($\frac{704 \mu g}{ml}$)	100 mM Vol. PBX
A	Spherotech 4.35 μm	10 mg	450 μl	322.6 μl	677.4 μl
B	Bangs 9803CN	2 mg	20 μl	283.0 μl	717 μl
C	Bangs 9500CN	2 mg	20 μl	199.0 μl	801 μl

- ① Add appropriate beads to a labeled 12x75 mm polypropylene tube.
- ② Wash bead 3x 1ml with 100 mM phosphate buffer pH 6.8 by adding 1ml buffer, vortexing and placing tube in Corning magnetic separator for 3 minutes.
- ③ Suspend beads in the volume of 100 mM phosphate buffer indicated above table.
- ④ Add the volume of CMV antigen specified in the table to the appropriate tube.
- ⑤ Place the capped tubes on an end-over-end rotator ON @ RT.
- ⑥ The next day place the tubes on a magnetic separator for 3 minutes. Pipet off & discard supernatant.
- ⑦ Wash 4x 1ml w/wash buffer by adding 1ml of wash buffer, vortexing and placing tubes in a Corning magnetic separator for 3 minutes.
- ⑧ In a similar manner wash 2x 1ml w/storage buffer.
- ⑨ Suspend the beads in 1ml of storage buffer and store @ 4°C.

Witnessed & Understood by me,

[Signature]

Date


Invented by

Recorded by

M. Watkins

Date

C



1840 Industrial Dr. Suite 270
Libertyville, Illinois 60048
Tel : (708) 680 8922
Fax : (708) 680 8927

TECHNICAL DATA

PRODUCT: SPHERO™ Carboxyl Magnetic Particles, 4.0-4.5 μm
(U. S. Patent No. 5,091,206)

CAT. NO.: CM-40-10

LOT NO.: 101

SIZE: 10 ml

PARTICLE CONC.: 2.5% w/v


PRESERVATIVE: 0.05% Sodium Azide*

STORAGE: Room Temperature

CAUTION: Do not freeze.

NOTE: To achieve optimum particle suspension, resuspend by vortexing before use.

SEM ANALYSIS: Magnification: 1000X. Mean Diameter: 4.35 μm



$$\text{Density} = 1.22 - 1.25 \text{ g/cc}$$

$$\% \text{ magnetite} = 12\%$$

$$\begin{aligned} \# \text{part./mL} &= \frac{6W \times 10^{12}}{\rho \pi \phi^3} \\ &= \frac{(6)(0.025)(10^{12})}{(1.235)(\pi)(4.35)^3} \\ &= 4.70 \times 10^8 \text{ part./mL} \end{aligned}$$

*W = grams/mL (soln)
 ϕ = diameter (μm)
 ρ = density (g/mL)*

***WARNING:** Sodium Azide can react with Cu and Pb in plumbing to form explosive metal azides. Flush this reagent down drains with copious amounts of water.

NOTE: FOR RESEARCH APPLICATIONS ONLY. NOT FOR DIAGNOSTIC USE.

Adsorption of HSV Antigen ^{To} 3.18 μ m Spherotech beads

Purpose: To adsorb Ross Southern Labs HSV antigen to 3.18 μ m magnetic beads from Spherotech.

Procedure

- ① Add 200 μ l (5mg, 2.5%) of Spherotech 3.18 μ m beads to a 12x75 mm polypropylene tube.
- ② Wash beads 3x 1 ml with 100 mM phosphate buffer pH 6.8 by adding 1 ml buffer, vortexing, placing tube in Corning magnetic separator for 3 minutes and pipetting off supernatant.
- ③ Suspend bead in 185 μ l of 100 mM phosphate buffer pH 6.8
- ④ Add 815 μ l of HSV antigen (Ross Southern Labs).
- ⑤ Cap the tube and place on a end-over-end rotator ON @ RT.
- ⑥ The next day place the tube on a magnetic separator for 3 minutes. Pipet off & discard supernatant.
- ⑦ Wash 4x 1 ml w/wash buffer by adding 1 ml of wash buffer, vortexing and placing tube in a Corning magnetic separator for 3 minutes.
- ⑧ In a similar manner, wash 2x 1 ml w/storage buffer
- ⑨ Suspend the beads in 1 ml of storage buffer and store at 4°C.



1840 Industrial Dr. Suite 270
Libertyville, Illinois 60048
Tel : (708) 690 8922
Fax : (708) 690 8927

TECHNICAL DATA

PRODUCT: SPHERO™ Carboxyl Magnetic Particles, 3.0-3.9 μm
(U. S. Patent No. 5,091,206)

CAT. NO.: CM-30-10

LOT NO.: 101

SIZE: 10 ml

PARTICLE CONC.: 2.5% w/v


PRESERVATIVE: 0.05% Sodium Azide*

STORAGE: Room Temperature

CAUTION: Do not freeze.

NOTE: To achieve optimum particle suspension, resuspend by vortexing before use.

SEM ANALYSIS: Magnification: 1000X. Mean Diameter: 3.18 μm



***WARNING:** Sodium Azide can react with Cu and Pb in plumbing to form explosive metal azides. Flush this reagent down drains with copious amounts of water.

NOTE: FOR RESEARCH APPLICATIONS ONLY. NOT FOR DIAGNOSTIC USE.

$$\begin{aligned} \text{Density} &= 1.22 - 1.25 \text{ g/cc} \\ \% \text{ Magnetite} &= .12\% \quad w = \text{grams/mL of sol'n} \\ \# \text{ part./mL} &= \frac{6W \times 10^{12}}{\rho \pi \phi^3} \quad \phi = \text{diameter } (\mu\text{m}) \\ &= \frac{(6)(0.025)(10^{12})}{(1.235)(\pi)(3.18)^3} \quad \rho = \text{density (g/mL)} \\ &= 1.20 \times 10^9 \text{ particles/mL} \\ \text{Surface Area} &= \frac{6\phi}{\phi \rho} = \frac{6\phi}{(3.18)(1.235)} = 15.3 \frac{\text{cm}^2}{\text{mg}} \end{aligned}$$

14

Adsorption of Rubella ^{To} 10 μ m Magnetic Sinter Beads

Project No. _____

Book No. _____

Date: _____

Purpose: To adsorb Rubella antigen at two different concentrations to 10 μ m magnetic sinter beads.

Procedure

Tube	Vol. Rubella Antigen	Vol. 100 mM PBS pH 6.8
------	----------------------	------------------------

A	104 μ L (5.2 μ g)	896 μ L
B	10.4 μ L (0.52 μ g)	989.6 μ L

- Wash 5 mg beads in tubes A & B, 3 x 1 ml of 100 mM phosphate buffer, pH 6.8 using magnetic separator (3 minutes).
- Suspend pellet in specified volume of phosphate buffer (see table).
- Add volume of rubella antigen specified in table.
- Place on end-over-end rotator O/N @ RT.
- Place on magnetic separator 3 minutes, discard ^{supernatant}.
- Wash 4 x 1 ml with wash buffer 3 min. of magnetic separation. ^{using}
- Wash 2 x 1 ml with storage buffer 3 min. of magnetic separation. ^{again using}
- Suspend in 1 ml of storage buffer.

To Page No. _____

Witnessed / Undertaken by: _____

Date: _____

Invented by: _____

Date: _____

(5)



Bio-Rad Labs
4000 Alfred Nobel Drive
Hercules, CA 94547
USA

SINTEF Applied Chemistry

Address:
N-7034 Trondheim,
NORWAY
Location:
Sem Sælands vei 2A
Telephone:
+47 73 59 28 73
Fax:
+47 73 59 69 95

Enterprise No.:
NO 948 007 029 MVA

Att: Dr. Mike Watkins

Your ref.:

Our ref.:

Direct line:
+4773592815

Trondheim,

MAGNETIC MICROSPHERES

Dear Dr. Watkins,

Please find enclosed 50mg of uncoated magnetic particles with the following specifications:

R-509: 10µm porous, superparamagnetic particles
surface area: 89m²/g
iron content: 17.9% Fe/g particles
(in the form of magnetite Fe₃O₄ and/or maghemite γ-Fe₂O₃)
magnetic susceptibility: 12·10⁻³ cgse

Density = 1.23 g/mL
particles/mL =

surface area (smooth) = $\frac{60}{(10)(1.23)} = 4.88 \frac{\text{cm}^2}{\text{mg}}$

We have several types of coated particles based on these uncoated beads, where the coating both serves as pore filler (→ compact, smooth surface) and as supplier for functional groups for ligand coupling. We can also design new coatings specially for your purpose. Shortly told, we can vary the surface area and the pore sizes, the surface chemistry, the Fe-content (→ the magnetic susceptibility) and the size.

Please use always our particle number R-509 in your further correspondence concerning these particles.

We are looking forward to hear about your experiences with these magnetic beads.

Yours sincerely
SINTEF Applied Chemistry

Ruth Schmid

Ruth Schmid
Research Scientist

cc

TITLE: Flow Multi (CMV + HSV + RUB) vs single (RUB) assay

Form Page No. _____

Purpose: To compare Rubella results in a single vs multiasay (HSV -, CMV, RUB) format.

Observations

- controls CN 6, CN 8, CN 12, CN 15 + 23 were tested with the Gull assay and found to have the following reactivities:

	<u>HSV</u>	<u>CMV</u>	<u>RUB</u>
CN 6	+	-	+
CN 8	+	-	+
CN 12	-	-	+
CN 15	-	-	+
23	-	+	-

the flow results are consistent with these reactivities.

- standard 5 gave a lower signal than standard 4 in both assay formats
- controls were lower than their reported value of 134.9, 14.4, 0.5 IU/ml for HI positive, low positive and negative controls, respectively.

Purpose: To use compare Rubella in single versus multitassay (i.e. HSV2, CMV and RUB).

Watkins

ASH

0282-28

Chemical

22

(01) 04-077

1

Methods

501

741

793 793

357

255

036 013

362 39

1

٧٠

**100 μ L sample (1/10 dilution with diluent)
add 100 μ L beads**

incubate 15 min. on vortex @ RT


**Add 750 μ L eluent, place on Corning magnetic separator 5 minutes
decant and let drain 1 minute on paper towels**

**Add 1000 μ L eluent, place on Corning magnetic separator 5 minutes
decant and let drain 1 minute on paper towels**

add 200 μ L anthrhuman IgG-PE(B) - staggered additions 3 minutes apart

incubate 15 min. on vortex @ RT

Read on Brite - staggered 3 minutes apart



		HSV			CMV			RUB		
Sample	Signal (Channels)	Rel. Linear Units	Total Counts	Signal (Channels)	Rel. Linear Units	Total Counts	Signal (Channels)	Rel. Linear Units	Total Counts	[Rub]
RUG 0							623	16	163	0.00
RUG 1							1080	129	186	7.97
RUG 2							1223	245	226	30.15
RUG 3							1304	352	115	95.41
RUG 4							1341	416	217	240.97
RUG 5							1226	248	172	31.19
# Pos							1259	268	167	46.79
o Pos							1100	141	188	9.41
Neg							704	24	266	0.21
CN 6							1182	204	209	19.65
CN 8							1174	196	189	18.18
Cn 12							1215	236	229	27.60
Cn 15							1313	367	199	114.44
23							875	51	187	1.51
RUG 0	374	5	741	470	8	652	613	16	258	0.00
RUG 1	501	10	635	479	9	447	1054	114	208	7.88
RUG 2	633	17	514	572	13	453	1189	220	309	30.55
RUG 3	741	28	572	698	23	378	1285	323	232	94.28
RUG 4	805	37	585	785	31	412	1330	396	244	242.36
RUG 5	793	35	524	931	66	406	1289	329	263	100.88
# Pos	634	17	443	1795	36	348	1241	265	180	49.94
o Pos	457	8	589	578	14	361	1063	130	211	10.11
Neg	363	5	484	700	23	442	702	24	296	0.22
CN 6	591	14	515	427	7	343	1122	155	213	14.31
CN 8	613	16	514	414	6	367	1156	181	267	19.69
Cn 12	359	5	535	425	7	365	1219	240	187	38.23
Cn 15	362	5	511	393	6	408	1289	329	178	100.88
23	346	5	516	844	45	348	836	43	246	1.21

ALLOWED CLAIMS (U.S. Patent Application No. 09/905,338)

21. (Amended) A composition comprising a plurality of solid-phase assay reagents selectively active in a plurality of assays each for a different analyte, each said solid-phase assay reagent comprising a binding species that is selectively active in a single assay and coupled to one of a plurality of microparticles of magnetically responsive material, the sizes of said microparticles varying in size over a range that is an aggregate of a plurality of subranges, each subrange distinguishable from other subranges of said aggregate by flow cytometry and by the binding species coupled thereto, said microparticles being suitable for use in a multiplex assay procedure that includes the use of flow cytometry.

22. A composition in accordance with claim 21 in which said range is a diameter from about 0.3 micrometers to about 100 micrometers.

23. A composition in accordance with claim 21 in which said range is a diameter of from about 0.5 micrometers to about 40 micrometers.

24.. A composition in accordance with claim 21 in which the standard deviation of the particle diameters of each said subrange is less than one third of the separation of the mean diameters of adjacent subranges.

25. A composition in accordance with claim 21 in which said microparticles have a porosity substantially less than macroporous.

26. A composition in accordance with claim 21 consisting essentially of from two to 100 binding species, each selectively active in a single assay relative to the remaining binding species.

27. A composition in accordance with claim 21 in which said microparticles are comprised of a combination of a polymer and a paramagnetic substance.

28. A composition in accordance with claim 27 in which said paramagnetic substance is a metal oxide.

29. A composition in accordance with claim 27 in which said polymer is formed from monomers including carboxylate groups to permit covalent bonding of assay binding members at the microparticle surface.

50. A composition according to claim 21 in which the microparticles are comprised of from about 1% to about 75% by weight of magnetically responsive material.

51. A composition according to claim 21 in which the microparticles are comprised of from about 2% to about 50% by weight of magnetically responsive material.

52. A composition according to claim 21 in which the microparticles are comprised of from about 3% to about 25% by weight of magnetically responsive material.

53. A composition according to claim 21 in which the microparticles are comprised of from about 5% to about 15% by weight of magnetically responsive material.

54. A composition according to claim 21 in which the microparticles are differentiable by size and by a differentiation parameter other than size.

55. A composition according to claim 54 in which the differentiation parameter other than size is one or more differentiation parameters selected from the group consisting of particle composition parameters, particle physical characteristics that affect light scattering, and dyes.

56. A composition according to claim 55 in which the differentiation parameter other than size is one or more differentiation parameters selected from fluorescence, colored dyes, light scatter, light emission, and absorbance.

57. A composition according to claim 54 in which the differentiation parameter other than size is one or more fluorescence parameters.

58. A composition according to claim 57 in which two or more fluorochromes are incorporated into each subrange of microparticles.